

# Lesson 6

## Instrument Control

### **You Will Learn:**


- **A. An overview of instrument control**
- **B. About GPIB communication and configuration**
- **C. About LabVIEW instrument drivers**
- **D. How to use instrument driver VIs**
- **E. About Virtual Instrument Software Architecture (VISA)**
- **F. How to use the VISA functions**
- **G. About serial port communication**
- **H. About waveform transfers**

# Instrument Control Overview



**Control any instrument if you know the following:**

- **Type of connector on the instrument**
- **Type of cables needed**
- **Electrical properties involved**
- **Communication protocols used**
- **Software drivers available**


**Instruments**




**Computer**



# GPIB Communication




# GPIB Hardware Specifications






- **Max cable length between devices = 4 m (2 m average)**
- **Max cable length = 20 m**
- **Max number of devices = 15 (2/3 powered on)**

# GPIB Software Architecture – Windows




# Configuring GPIB Board and Instruments

- **Measurement & Automation Explorer (MAX)**







# Configuring GPIB Board and Instruments






# Instrument Drivers

- More than 2,200 LabVIEW instrument drivers from more than 150 vendors
- Programming simplified to high-level Application programming interface (API)




# Installing and Finding Instrument Drivers

- NI Web site or Instrument Driver CD
- Install the instrument driver VI Library into LabVIEW\instr.lib directory
- Access drivers from Instrument I/O » Instrument Drivers subpalette




# Instrument Drivers

- Instrument drivers have a similar hierarchy




**Instrument Driver VI Tree**

# Instrument Driver Inputs and Outputs



- **VISA Sessions**
  - A connection or link to a specific instrument
  - Created after instrument is initialized
  - Used throughout VI whenever you communicate with that specific instrument
- **Error Clusters**


# VISA




## VISA VIs and Functions

The lower layer of functions in the LabVIEW instrument driver VIs that communicates with the driver software uses the Virtual Instrument Software Architecture (VISA). The VISA architecture enables instrument drivers to talk to different hardware interfaces using one standard application programming interface. You can use LabVIEW VISA functions and VIs to communicate with an instrument that does not have a driver or that has a driver that implements only a subset of the total capabilities of the instrument. This lesson shows how you can use the VISA Write and VISA Read functions to communicate with GPIB instruments.


# VISA Resource Name



- **Exact name and location of the instrument**
- **Use the VISA Resource Name control (like the DAQ Channel Name control)**
- **You can specify the full resource name or the VISA Alias**



# Putting It All Together




- Initialize instrument
- Do operation(s)
- Close instrument
- Check error status



# VISA

Interface	Syntax
Asynchronous	ASRL[ <i>board</i> ] [::INSTR]
GPIB	AGIP[ <i>board</i> ]:: <i>primary address</i> :: <i>secondary address</i> [::INSTR]
VXI instrument through embedded or MXIbus controller	VXI[ <i>board</i> ]::VXI <i>logical address</i> [::INSTR]
GPIB-VXI controller	GPIB-VXI[ <i>board</i> ] [::GPIB-VXI <i>primary address</i> ]::VXI <i>logical address</i> [::INSTR]



syntax for the VISA resource name varies depending on the interface to the instrument. The table on this slide identifies the syntax for the VISA resource name for the asynchronous, GPIB, VXI, and GPIB-VXI interfaces. You can use Measurement & Automation Explorer (MAX) to determine what resources and instrument addresses are available. Refer to Exercise 9-1 in this program for more information about how to use MAX.


When you initialize the communication channel to an instrument, you must know the VISA resource name or instrument descriptor. A resource is an instrument or interface and the instrument descriptor is the exact name and location of a resource in the following format:

*Interface Type*[*board index*]::*Address*::INSTR

For example, GPIB::2::INSTR is the instrument descriptor for a GPIB instrument at address 2. Use the VISA Resource Name control (shown above) located on the I/O subpalette of the **Controls** palette to supply the VISA resource name from the front panel controls of a VI.



# Virtual Instrument Software Architecture



- **Standard API for programming instruments**
  - **Platform independent**
  - **Interface independent**
  - **Interface to Serial, GPIB, and VXI instruments**
  - **Uses VISA.DLL**

# VISA Terminology

- Resource—Instrument, Serial Port, or Parallel Port
- Session—Connection to a Resource
- Instrument Descriptor—Resource location

- **Format: Interface**  
**Type::Address::INSTR**

- **Examples:**


GPIB0:1::INSTR
GPIB0:4::INSTR
GPIB0:10::INSTR
ASRL1::INSTR
ASRL2::INSTR
ASRL3::INSTR
ASRL10::INSTR

# Instrument Descriptor Syntax


- Resource Name contains interface info
- VISA Aliases also work

Interface	Resource Name Grammar
Serial	ASRL[board][::INSTR]
GPIB	GPIB[board]:: <i>primary address</i> [::INSTR]
VXI	VXI[board]:: <i>VXI logical address</i> [::INSTR]
GPIB-VXI	GPIB-VXI[board][:: <i>GPIB-VXI primary address</i> ]:: <i>VXI logical address</i> [::INSTR]


# VISA Functions




## VISA Write




## VISA Read



# Serial Communication




- Popular means of communication between computer and peripheral device
- Data sent one bit at a time across the cable
- Used for low transfer rates or long distances
- Only a cable is needed since most computers have at least one available serial port



## Serial Port Communication

Use the VIs and functions located on the **Serial** subpalette for serial port communication. You used the VISA Write and VISA Read functions on this palette for GPIB communication in Exercise 9-4. The VISA Read and VISA Write functions work with any type of instrument communication, including serial. However, because serial communication requires you to configure extra parameters, you must start the serial port communication with the VISA Configure Serial Port VI. This lesson shows how you can use the VISA Configure Serial Port VI to communicate with an instrument through a serial interface.

# Serial Communication




## Terminology

- **Baud rate – bits per second**
- **Data bits – inverted logic and LSB first**
- **Parity – optional error-checking bit**
- **Stop bits – 1, 1.5, or 2 inverted bits at data end**
- **Flow control – hardware and software handshaking options**

# Serial Hardware Connection


- **RS-232**
  - Most PCs
  - DCE or DTE configurations
  - 8-pin or 25-pin
  - Single-ended
- **RS-422**
  - Macintosh
  - 8-pin
  - Differential
- **RS-485**
  - Multidrop
  - Industrial Automation




Pin		DTE	DCE
1	DCD	Input	Output
2	RxD	I	O
3	TxD	O	I
4	DTR	O	I
5	Com	-	-
6	DSR	I	O
7	RTS	O	I
8	CTS	I	O
9	RI	I	O



# Serial VIs and Functions




- Found in Serial subpalette under Instrument I/O
- Based on VISA functions
- Serial VIs and functions also work with parallel port communication



**VISA Configure Serial Port**

# Serial I/O Example



- Initialize the serial port settings
- Write commands to the device
- Read device response
- Check for errors


# Serial I/O Example

National Instruments stepper motor

Front Panel

Abort deceleration 500	Acceleration 150	VISA out .
Target position 1000000000	starting position -1000000000	Record serial data (0=off;1=on) 0
Abort action code 0=disable drive 1=stop and servo 2=decelerate and stop 2	Setup Cluster VISA in . serial port ----- Address node 0 checksum NO baud rate 9600 log serial data NO serial timeout 0 Wait_after_send 0-255 [ mSec] 0 max.re-sends 0	error OUT status code ✓ 0 source
Speed 4000		
Deceleration 150		
error in (no error) status code ✓ 0 source		

## Block Diagram



# Summary

- **LabVIEW can communicate with any instrument that connects to your computer if you know the interface type**
- **Use the Measurement & Automation Explorer (MAX) to detect, configure, and test your GPIB interface and instruments**
- **An instrument driver eliminates the need for your to have detailed knowledge of the specific strings used by an instrument**
- **Instrument Library – more than 650 instruments supported**
- **Instrument driver VIs share a common hierarchy and come with an example to help you get started**
- **VISA a standard protocol for using multiple types of I/O and instrument driver development**
- **Serial library contains functions for serial communication**
- **You need to know the format of the returned data string in order to convert it to the correct values**